# ProteinPurification

*NellySelem*

*May 18, 2017*

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Flask each one with a 100ml were prepared. columns were washed with 20ml of this tubes at 2000rpm 5 min  
Imidazol concentration=x

v: (x=)100ml/500mM

Flask	BufferB (mM)	BufferB (ml)	Buffer A ml	Buffer B
1	50	10	90	10
2	100	20	80	20
3	150	30	70	30
4	300	60	40	60
5	500	100	0	100

Flask must be kept cold.

7. Recover the fractions that contains the protein that elutes approximately at 150-200 mM imidazol. (TrpC, TrpD) or 300mM imidazol (must PriA). Collect the relevant fractions. At this point a protein gel can be run to check protein elution.
8. Concentrate the protein using AMICON tubes 10kDa (all molecules below this weight will be discarded, molecules with weight above 10kDa such as PriA will be kept inside the tube) until an approximate concentration of  $90 \frac{\mu g}{\mu l}$  (1ml vol approx). Because imidazol > 200mM reacts with Bradford, protein on amicon must be washed with 10 ml of buffer A to eliminate imidazol (4ml each time on amicon tube). Finally, after protein concentration, 100 ul sample must be recovered with 900 ul Buffer A (free of imidazol) (8000 rpm 10 min) when AMICON are stored for long periods leave on ethanol at 20%, to reuse them throw ethanol out and wash many times with milliQ water.
9. Check purity by SDS-PAGE gel  
on *S. coelicolor* PriA molecular weight is around 25.7 kDa (length 240 aminoacids), on *T. maritima* TrpC molecular weight is around 27.5 kDa (length 242 aminoacids)
10. Add glycerol at a final concentration of 15% and quantify using Bradford. Store at -80C

*TrpC important before step 5*

After centrifuge step 4 :

Heat supernatant by 15 min at 65C

Centrifuge at 10000rpm by 20min. Recover the supernatant.

Continue with steps 5-10 as in trpD purification.

On protein activity is more important to be pure than to have protein abundance.